

Does Vessel Prep Make a Difference?

AVAFLEX Study 6 Month Results

JOHN ROSS, MD

Medical Director, The Regional Medical Center

Orangeburg, SC

Resistant Stenosis & Restenosis in AVF/AVG


Problem:

- Failure of fistulas remains a significant issue and costs the healthcare system \$2B annually¹

Results:

(Tretotola et al, JVIR, 2005):

- 55% of Lesions require pressures greater than 15 atm to efface the waist
- Primary patency AVF/AVG through 12 months: 48-80% 24 months: 50-80% (Holden CIRSE 2019)




¹ USRDS: USRDS 2010 Annual Data Report: Atlas of Chronic Kidney Disease and End-Stage Renal Disease in the United States, Bethesda, MD, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases, 2010

Primary Patency of PTA in AV Access


Study Year	Primary Patency Reported Rates			Source
	6 months	12 months	24 months	
2000	67% Forearm AVF	51% Forearm AVF	37% Forearm AVF	Turmel-Rodrigues et al. 2000, n=220 dysfunctional AVF
2004	61% After PTA in AVFs	42% After PTA in AVFs	35% After PTA in AVFs	Bountouris et al., 2014, n=159 PTA AVF's

Most other studies have reported similar primary patency rates of ~ 40–50 % at 1 year¹⁻³



1. Rajan DK, Bunston S. Dysfunctional autogenous hemodialysis fistulas: outcomes after angioplasty – are there clinical predictors of patency? Radiology. 2004;232:508–15.
2. Manninen HI, Kaukanen ET. Brachial arterial access: endovascular treatment of failing Brescia-Cimino hemodialysis fistulas – initial success and long-term results. Radiology. 2001;218:711–8.
3. Heye S, Maleux G. Factors influencing technical success and outcome of percutaneous balloon angioplasty in de novo native hemodialysis arteriovenous fistulas. Eur J Radiol. 2012;81:2298–303. doi: 10.1016/j.ejrad.2011.09.004.

Lutonix DCB in AVF vs PTA



Conclusions Paclitaxel-coated balloon angioplasty did not meet the primary effectiveness endpoint at 180 days compared to conventional angioplasty. Both arms demonstrated equivalent safety.


Scott Trerotola, Jeffrey Lawson, Prabir Roy-Chaudhury, and Theodore Saad, for the Lutonix AV Clinical Trial Investigators. **Randomized Trial of Drug Coated Balloon Angioplasty in Failing AV Fistulas.** CJASN doi: 10.2215/CJN.14231217

Why Vessel Prep Prior to Angioplasty in AV Access?

- PTA has significant limitations:
 - Elastic recoil
 - Resistant hyperplasia – even with ultra high-pressure balloons
- Vessel prep along the entire length of stenosis:
 - Improves vessel compliance
 - Eliminates the need for high pressure balloons
 - Minimizes barotrauma that may extend primary patency or time to reintervention



FLEX Vessel Prep™ System – AV



FLEX Vessel Prep™ System - Dual MOA



Longitudinal Micro-incisions Combined With Radial Expansion Force

Controlled, parallel micro-incisions

Radial expansion force

3 protective “atherotome-mounted” skids


- Proprietary depth control “flexes” to guide along the contour of the lesion



DCB Studies In AVF/AVG

Lutonix IDE Study:

- Multicenter
- Randomized
- Prospective




IN.PACT AV Access IDE Study:

- Multicenter
- Randomized
- Prospective




AVAFLEX Study:

- Single Center
- Single Arm
- Prospective



Lutonix AV IDE - 6 Month Data

Target Lesion Primary Patency at 6 Months (210 Days)



**6 Month Lutonix
Target Lesion
Patency: 71%**

Source: lutonixdcb.com/av-treatment

Trerotola SO, Lawson J, Roy-Chaudhury P, Saad TF, et al.

Drug Coated Balloon
Angioplasty in Failing AV
Fistulas: A Randomized
Controlled Trial. ClinJ Am
SocNephrol2018;13:1215-1224.

AVAFLEX Study

FLEX - Drug Coated Balloon Dialysis ACCESS Stenosis Study

Start Date: May 2018

Lead Investigator: John Ross, MD, FACS

Investigator Location: Dialysis Access Institute
Orangeburg, South Carolina

Study Objective: To evaluate the efficacy of the FLEX Vessel Prep™ System with a Drug Coated Balloon (DCB) in the maintenance of arteriovenous access

Similar Key Inclusion/Exclusion Criteria to Lutonix IDE Trial.



AVAFLEX Demographics


Characteristic	Mean \pm SD or n (%) (N=34)
Gender	
Female	13 (38.2%)
Male	21 (61.8%)
Age	60.4 \pm 11.0
Smoker	3 (8.8%)
Hypertension	33 (97.1%)
CAD	13 (38.2%)
Congestive HF	14 (41.2%)
Diabetes	17 (50.0%)

Characteristic	Mean \pm SD or n (%) (N=34)
Type	
AVF	31 (91.2%)
AVG	3 (8.8%)
# Prior Intervention	
0	12 (35.3%)
1	6 (17.6%)
2	8 (23.5%)
3	2 (5.9%)
4	1 (2.9%)
5	1 (2.9%)
≥ 6	4 (11.8%)


AVAFLEX Procedural Characteristics

Procedural Characteristics	
Target Vessel Diameter (mm)	8.7 ± 1.1
Target Lesion Length (cm)	2.3 ± 2.6
Cephalic Arch Lesion	52.9% (18/34)
Pre-Stenosis	79.9 ± 9.9
Post FLEX Lumen Gain	19.1 ± 20.1
Post Procedure Stenosis	7.4 ± 6.6
Device Success*	100% (34/34)

*Device Success is defined as the successful delivery to the target lesion, deployment of the treatment element, and retrieval at the index procedure.




Pre-Fistulagram



Post-FLEX + Lutonix DCB Fistulagram

6 Month Target Lesion Patency




AVAFLEX
(FLEX + Lutonix DCB)
Target Lesion Patency:
84.1%

Lutonix DCB
Target Lesion Patency:
71%

Lutonix POBA
Target Lesion Patency:
63%

AVAFLEX Data

Re-Weighted per Lutonix Lesion Population




Time to target vessel re-intervention by device – FLEX re-weighted per Lutonix Lesion Location Population
Cephalic arch strata weighted 18.7% (versus original population percentage of 52.9%)
Other target vessels weighted 81.3% (versus original population percentage of 47.1%)

Summary

AVAFLEX study demonstrates vessel prep with the FLEX Vessel Prep™ System improves 6 Month Target Lesion Patency prior to drug coated balloon or PTA.

AVAFLEX 6-Month Target Lesion Primary Patency is 84.1% (vs. 71% Lutonix AV IDE) with no procedural issues

In this challenging patient population, vessel prep has the potential to improve long-term patency due to:


- Reduced elastic recoil
 - Improved ability to dilate in resistant hyperplastic stenoses
 - Enhanced drug-uptake through micro-incisions
- 

FLEX-AV Registry

Objective: To document the endovascular intervention approaches and outcomes when the FLEX Vessel Prep™ System is utilized in a hospital per the institution's standard practice and at 6, 9, and 12 months following treatment.

Hypothesis: The FLEX VP System combined with Balloon Angioplasty will result in Anatomic Success of the index procedure and extend the period of target lesion primary patency as compared to balloon angioplasty alone as demonstrated with historic controls.

- **15-20 Registry Sites**
- **Western Institutional Review Board (WIRB) Approved**
- **Study initiation Start Date: September 2019**



John Aruny MD.

**Principal
Investigator**

Dialysis Access Institute
Orangeburg, SC

Thank You!