



# FLEX AV Registry 12-month Results





# FLEX-AV Registry

- **Study Design:** Multi-Center, Single-Arm, Prospective Study
- **Study Population:** Patients with AVF/AVG stenosis, eligible for FLEX/Angioplasty treatment.
- **Enrollment:** 114
- **Sites:** 8 Sites
- **Follow-Up:** Phone call to the patient and/or dialysis center at 6, 9, and 12 months.

# FLEX-AV Registry: Enrollment Overview

- 114 patients treated
- Total number of lesions treated, n= 148
- 8 clinical sites



# Patient Demographics & Characteristics: Real-World, No “Cherry-picking”

Variable	114 Subjects
Age (years)	63.3 ± 12.7 (114) 31.0-88.0
Gender	
Female	61/114 (53.5%)
Male	53/114 (46.5%)
Race	
American Indian or Alaska Native	2/114 (1.8%)
Asian	1/114 (0.9%)
Black or African American	75/114 (65.8%)
White	36/114 (31.6%)
Smoking History	
Current	17/114 (14.9%)
Never	60/114 (52.6%)
Past	37/114 (32.5%)

114 patients, 8 sites, All-comers Registry	
Prior AV Access Interventions	4.9 ± 5.8; highest number: 29
Cephalic Arch %	22%
Target Lesion Length (mm)	21 ± 25
Maximum PTA Balloon Pressure needed (atms)	15.2 ± 5.9
<i>No limitations on Lesion Length</i>	
<i>No limitations on Prior Interventions</i>	




## **Functional Primary Patency 12 months**

**Freedom from Target Lesion Reintervention**


# 12-month Target Lesion Primary Patency Analysis

- 112/114 target lesions were considered for the Kaplan-Meier estimate
  - 2 subjects did not contribute to any follow-up analysis
- Patients analyzed in following cohorts
  - AVF with PTA alone
  - AVG with PTA alone
  - Cephalic Arch target lesions with PTA alone
- Statistical Analysis
  - ✓ Functional Patency estimation via Kaplan-Meier analysis at the close of the 6-month FU visit
  - ✓ Freedom from TLR (time) was determined using restricted mean survival time analysis (RMST) restricted to 270 days

# Patient Cohorts





# 12-month Target Lesion Primary Patency FLEX + PTA



\* Kaplan Meier estimates are only provided when there are patients remaining under observation at 372days. The 9-month TLPP = 41.2%.



# 12-month Target Lesion Primary Patency FLEX+PTA in AVF



# 12-month TLPP Literature Comparisons

## PTA in AVF Results

Published Results	FLEX Registry FLEX + PTA	Liao, et al <sup>1</sup>	Rajan, et al <sup>2</sup>	Ng, et al <sup>3</sup>	Hu, et al <sup>4</sup>
12-month TLPP AVF	49.0% (n=34)	31.5% (n= 273)*	26% (n = 53)	0% - 21.2% (n=143)*	47.2% (n=341)*

\*restricted to AVF studies only



<sup>1</sup>Liao M-T, Chen M-K, Hsieh M-Y, Yeh N-L, Chien K-L, Lin C-C, et al. Drug-coated balloon versus conventional balloon angioplasty of hemodialysis arteriovenous fistula or graft: A systematic review and meta-analysis of randomized controlled trials. PLOS One; 2020 15(4).

<sup>2</sup>Rajan D., et al., Dysfunctional Autogenous Hemodialysis Fistulas: Outcomes after Angioplasty – Are There Clinical Predictors of Patency? Radiology. Sept 2004.

<sup>3</sup> Ng B, Fugger M, Onakpoya IJ, et al. Covered stents versus balloon angioplasty for failure of arteriovenous access: a systematic review and meta-analysis. BMJ Open 2021;11

<sup>4</sup> Hu H, Tan Q, Wang J, Liu Y, Yang Y, Zhao J. Drug-coated balloon angioplasty for failing haemodialysis access: meta-analysis of randomized clinical trials. Br J Surg. 2021 Nov 11;108(11):1293-1303.

# 12-month Target Lesion Primary Patency FLEX + PTA in Cephalic Arch



# 12-month TLPP Literature Comparisons PTA in Cephalic Arch Results

Published Results	FLEX Registry FLEX + PTA	D’Cruz, et al. <sup>1,5</sup>	Tng et al. <sup>2</sup>	Vasanthamohanm, et al. <sup>3</sup>	Miller, et al. <sup>4</sup>
12-month Functional Patency Cephalic Arch	<b>59.7% (n=14)</b>	9.5%* (n=146)  *Pooled multi-studies (0-39.6%)	33.9% (n=59)	0%-23% (n= 13-24)  * multi-studies small sample sizes	11% (N=50)*  * Historical controls

<sup>1</sup> D’Cruz RT, Leong SW, Syn N, et al. Endovascular treatment of cephalic arch stenosis in brachiocephalic arteriovenous fistulas: a systematic review and meta-analysis. J Vasc Access 2019; 20: 345.



<sup>2</sup>Tng RK, et al., Treatment of cephalic arch stenosis in dysfunctional arteriovenous fistulas with paclitaxel-coated versus conventional balloon angioplasty, CVIR Endovascular, (2021) f:80.

<sup>3</sup> Vasanthamohanm, L., et al. The Management of Cephalic Arch Stenosis in Arteriovenous Fistulas for Hemodialysis: A Systematic Review Cardiovasc Intervent Radiol (2015) 38:1179–1185

<sup>4</sup> Miller GA, Preddie DC, Savransky Y, Spergel LM. Use of the Viabahn stent graft for the treatment of recurrent cephalic arch stenosis in hemodialysis accesses. J Vasc Surg. 2018.

<sup>5</sup> Beathard et al., End Points for Interventional Studies for AV Access, Clin J Am Soc Nephrol 13: 501–512, March 2018.

# 12-Month Functional Patency FLEX + PTA in AVG



\*Kaplan-Meier estimates are only provided when there are patients remaining under observation at 372

# 12-month TLPP Literature Comparisons PTA in AVG Results

Published Results for PTA of AVGs	FLEX Registry FLEX + PTA	Yang, et al <sup>1</sup>	Liao, et al. <sup>2</sup>	Ng, et al. <sup>3</sup>
12-month Functional Patency AVG	--%* (n=25)	7.8% (N=49)	9% (n=22)	0% - 24.8% (n=339) <small>*Meta-analysis of RCTs of covered stents – AVGs</small>

**9 Month  
41.2% (n=18)**

\*Kaplan-Meier estimates are only provided when there are patients remaining under observation at 372

<sup>1</sup> Yang HT, Yu SY, Su TW, Kao TC, Hsieh HC, Ko PJ. A prospective randomized study of stent graft placement after balloon angioplasty versus balloon angioplasty alone for the treatment of hemodialysis patients with prosthetic graft outflow stenosis. *J Vasc Surg.* 2018 Aug;68(2):546-553.

<sup>2</sup> Liao M-T, Chen M-K, Hsieh M-Y, Yeh N-L, Chien K-L, Lin C-C, et al. Drug-coated balloon versus conventional balloon angioplasty of hemodialysis arteriovenous fistula or graft: A systematic review and meta-analysis of randomized controlled trials. *PLOS One;* 2020 15(4).

<sup>3</sup>Ng B, Fugger M, Onakpoya IJ, et al. Covered stents versus balloon angioplasty for failure of arteriovenous access: a systematic review and meta-analysis. *BMJ Open* 2021;11


# Procedure Complications Reported

- **No Serious Adverse Events**
- 4.3% Procedure complications reported during the study

Classification	Reported Complication	Number	Reported Cause	Treatment	FLEX %	JVIR Quality Improvement Guideline <sup>1</sup> Thresholds
Major	None	0	N/A	N/A	0%	2% (AVF) 7% (AVG)
Minor	Dissection, Grade B,C	4	Angioplasty	1 – No treatment 3 – Balloon inflation	4.3%	8% (AVF) 4% (AVG)
	Balloon burst	1	Balloon Rupture	Embolectomy		

<sup>1</sup> Daruiushnia, S. et al. Quality Improvement Guidelines for Percutaneous Image-Guided Management of the Thrombosed or Dysfunctional Dialysis Circuit, J VascInterv Radiol 2016;27:1518–1530

# Access Circuit Primary Patency at 12 months



\*Kaplan-Meier estimates are only provided when there are patients remaining under observation at 372



# Comparison ACPP

Published Results for ACPP	FLEX Registry FLEX + PTA	Dolmatch 2023 <sup>1</sup>	Holden 2022 <sup>2</sup>	Fong 2021 <sup>3</sup>	Haskel 2016 <sup>4</sup>
12-month ACPP in AVF	<b>43%</b>	17.7% (n=138)	32.4% (n=160)	29.8% (n=424)	n/a
12-month ACPP in AVG	<b>---%*</b>	n/a	n/a		11% (n=132)

# Conclusions

- The FLEX AV Registry 12-month outcomes demonstrate sustained patency across most subjects and impressive results specifically in the Cephalic Arch.
- Result highlights for FLEX + PTA:
  - 49% patency at 12 months
  - 60% patency at 12 months
  - 41.2% patency for all AVG patients at 9 months
  - AVG Days, Freedom from TLR
    - All (AVF & AVGs) Flex + PTA **256.6 (8.5 months)**
    - All AVF, FLEX + PTA **267.3 (8.9 months)**
    - All CA + PTA **267.7 (8.9 months)**
    - All AVG, Flex + PTA **228.5 (7.6 months)**
  - Lower PTA max inflation pressures 15.2 ATM
  - No observed SAEs